
Rabbit Anti-Monkey IgG (H&L) - HRP

Catalog: BGT-ANS-07010 Host: Rabbit

Description:

Rabbit Polyclonal Secondary Antibody to Monkey IgG (H&L) HRP

Product:

1 mg/ml, liquid in 0.01M Phosphate Buffered Saline, pH 7.2, containing 1% BSA, 50% glycerol, 0.02% Proclin

Preservative:

0.02% Proclin/Sodium Azide

Classification:

Secondary Antibodies

Purification&Purity:

The antibody was isolated from antisera by immunoaffinity chromatography using antigens coupled to agarose beads.

Applications:

E (1/5000 - 1/20000), WB (1/5000 - 1/20000), IH (1/500 - 1/1000), IC (1/500 - 1/1000)

Storage&Stability:

Shipped at 4 °C. Upon delivery aliquot and store at -20 °C for one year. Avoid freeze/thaw cycles.

Specificity:

By immunoelectrophoresis and ELISA this antibody reacts specifically with Monkey IgG. No antibody was detected against non immunoglobulin serum proteins.

Expiration Date:

Six months from date of reconstitution.

Note:

For research use only, not for use in diagnostic procedure.
