

Theophylline Antibody

Summary

| | |
|-------------------------------------|--|
| Catalog No. | BGT-ANT-03036 |
| Host species | Mouse Source: Ascites |
| Predicted Species Reactivity | Theophylline |
| Clone | B873M |
| Form | Purified, Liquid |
| Storage buffer | PBS, pH 7.0 |
| Preservative | 0.02% Sodium azide |
| Purity | >95% pure (SDS-PAGE) |
| Description of Target | Monoclonal Antibody to Theophylline |
| Clonality | Monoclonal |
| Isotype | IgG1 |
| Immunogen | Theophylline-8-BSA |
| Specificity | Theophylline |
| Purification | Purified |
| Reconstitution and Storage | Short term store at 2 to 8°C. Long term store at -20°C. Avoid multiple freeze/thaw cycles. |
| Warnings | This product contains sodium azide, which has been classified as Xn(Harmful), in European Directive 67/548/EEC in the concentration range of 0.1 - 1.0 %. When disposing of this reagent through lead or copper plumbing, flush with copious volumes |