

Datasheet

Phenytoin polyclonal antibody

Catalog Number: BGT-ANT-03041

Regulatory Status: For research use only (RUO)

Product Description: Sheep polyclonal antibody raised against Phenytoin.

Immunogen: Phenytoin conjugated with BTG.

Host: Sheep

Applications: ELISA

Form: Liquid

Recommend Usage: ELISA (1.25 ug/mL)

The optimal working dilution should be determined by the end user.

Storage Buffer: In 20 mM phosphate buffer, 150 mM NaCl, pH 7.2 (0.09% sodium azide)

Storage Instruction: Store at -20°C.

Aliquot to avoid repeated freezing and thawing.