

---

## PRODUCT SPECIFICATIONS

Product Name	Feline Serum amyloid A antibody
Catalog#	BGT-ANT-07068
Source	Mouse
Format	Monoclonal antibody
Formulation	Supplied as a 0.22 $\mu$ m filtered solution of PBS, pH7.4
Specificity	Feline Serum amyloid A (SAA) protein
Isotype	IgG1
Concentration	Lot specific( $\pm$ 20%)
Purification	Protein A or G affinity purification chromatography
Purity	$\geq$ 90% as determined by SDS-PAGE
Appearance	Transparent clear
Application	Elisa, Lateral flow, Chemiluminescence or other applications
Preservative	0.05% Proclin 300
PI	NA
Shipping	The product is shipped on dry ice pack. Upon receipt, store it immediately at -20°C or lower.
Storage	Short term (up to 7 days) store at 2-8°C. Long term, store at -20°C. Avoid multiple freeze/thaw cycles.
Expiration	Store at -20°C for three years