



Product name: PerCP Mouse IgG1, κ Isotype Control

Catalog: BGT-ANT-07091

Form: Liquid

Isotype: Mouse IgG1, κ

Size: 25 μ g/100 μ g

Host: Mouse

Concentration: 0.2 mg/mL

Conjugation: PerCP

Buffer: Phosphate buffered solution, pH 7.2, containing 0.09% sodium azide and 1% BSA

Recommended Use:

Each lot of this antibody is quality control tested by flow cytometric analysis as negative control. Use at concentrations comparable to those of the specific antibody of interest.

Storage:

Keep as concentrated solution.

Store at 2~8°C and protected from prolonged exposure to light. Do not freeze.

This product is guaranteed up to one year from purchase.