

Datasheet

WDR42A monoclonal antibody

Catalog Number: BGT-ANT-23684

Regulatory Status: For research use only (RUO)

Product Description: Mouse monoclonal antibody raised against a full-length recombinant WDR42A.

Clone Name: 4G8

Immunogen: WDR42A (AAH13107, 1 a.a. ~ 273 a.a) full-length recombinant protein with GST tag. MW of the GST tag alone is 26 KDa.

Sequence:

MSSKGSSTDGRTDLANGSLSSSPEEMSGAEEGRETS
SGIEVEASDLSLSLTGDDGGPNRTSTESRGTETESSG
EDKSDSMEDTGHYSINDENRVHDRSEEEEEEEEEEE
EEQPRRRVQRKRANRDQDSSDDERALEDWVSSETSA
LPRPRWQALPALRERELGSSARFVYEACGARVQVQR
FRLQHGLEGHTGCVNTLHFNQRGTWLASGDDLKVV
VWDWVRRQPVLDVFESGHKSNVFQVRQGSIIATERIRH
ELEKSELQHGLGADSELL

Host: Mouse

Interspecies Antigen Sequence: Mouse (95); Rat (96)

Reactivity: Human

Applications: ELISA, IF, S-ELISA, WB

Isotype: IgG2a Kappa

Storage Buffer: In 1x PBS, pH 7.4

Storage Instruction: Store at -20°C or lower. Aliquot to avoid repeated freezing and thawing.

Entrez GeneID: 50717

Gene Symbol: WDR42A

Gene Alias: DKFZp781G1096, FLJ35857, H326, MGC117276, MGC118891, MGC99640

Gene Summary: This gene encodes a WD repeat-containing protein that interacts with the Cul4-Ddb1 E3 ligase macromolecular complex. Multiple alternatively spliced transcript variants have been found for this gene. [provided by RefSeq]