

# CD62P/SELP Monoclonal Antibody

## Summary

---

<b>Catalog No.</b>	BGT-ANT-35055
<b>Alternative Names</b>	Leukocyte-endothelial cell adhesion molecule 3, Platelet activation dependent granule-external membrane protein, Granule membrane protein 140, GMP-140, CD62P, SELP, GRMP, PADGEM, CD62 antigen-like family member P, LECAM3, P-selectin, GMRP
<b>Host species</b>	Human
<b>Conjugation</b>	Unconjugated
<b>Species reactivity</b>	Human
<b>Form</b>	Liquid
<b>Storage buffer</b>	0.01M PBS, pH 7.4.
<b>Concentration</b>	6.06 mg/ml
<b>Purity</b>	>95%
<b>Clonality</b>	Monoclonal
<b>Isotype</b>	IgG2, kappa
<b>Applications</b>	ELISA, FCM
<b>Purification</b>	Protein A or G purified.
<b>Endotoxin level</b>	Please contact with the lab for this information.
<b>Accession</b>	P16109
<b>Stability and Storage</b>	Use a manual defrost freezer and avoid repeated freeze-thaw cycles. Store at +4°C short term (1-2 weeks). Store at -20 °C 12 months. Store at -80°C long term.