


Ganglioside GM1 Monoclonal Antibody

Summary

Catalog No.	BGT-ANT-35056
Alternative Names	Ganglioside GM1, Ganglioside A2, Ganglioside G4, Ganglioside GGtet1, Ganglioside GI, Monosialoganglioside GM1
Host species	Mouse
Form	Liquid
Storage buffer	0.01M PBS, pH 7.4.
Concentration	1 mg/ml
Purity	>95% by SDS-PAGE
Clonality	Monoclonal
Isotype	IgG2a, kappa
Applications	ELISA
Target	Ganglioside GM1
Purification	Protein A or G purified.
Endotoxin level	Please contact with the lab for this information.
Stability and Storage	Use a manual defrost freezer and avoid repeated freeze-thaw cycles. Store at +4°C short term (1-2 weeks). Store at -20 °C 12 months. Store at -80°C long term.
Note	For research use only.

Data Image



The purity of this product is >95% as determined by SEC-MALS.

SEC-HPLC