

Clostridium botulinum botA/BoNT Antibody

Summary

Catalog No.	BGT-ANT-35235
Clone ID	F1-40
Host species	Mouse
Species reactivity	Clostridium botulinum (strain Hall / ATCC 3502 / NCTC 13319 / Type A)
Form	Liquid
Storage buffer	0.01M PBS, pH 7.4.
Concentration	1 mg/ml
Purity	>95% as determined by SDS-PAGE.
Clonality	Monoclonal
Isotype	IgG2a
Applications	ELISA, WB
Target	botA, Botulinum neurotoxin type A, BoNT/A, Bontoxilysin-A, BOTOX
Purification	Protein A/G purified from cell culture supernatant.
Expression system	Mammalian Cells
Accession	P0DPI1
Stability and Storage	Use a manual defrost freezer and avoid repeated freeze-thaw cycles. Store at 4°C short term (1-2 weeks). Store at -20°C 12 months. Store at -80°C long term.