


dsDNA Antibody

Summary


Catalog No.	BGT-ANT-35343
Alternative Names	DNA, dsDNA
Host species	Mouse
Form	Liquid
Storage buffer	0.01M PBS, pH 7.4.
Concentration	0.68 mg/ml
Purity	>95%
Clonality	Monoclonal
Isotype	IgG2a, kappa
Applications	Dot, ELISA
Purification	Protein A or G purified from cell culture supernatant.
Endotoxin level	Please contact with the lab for this information.
Stability and Storage	Use a manual defrost freezer and avoid repeated freeze-thaw cycles. Store at +4°C short term (1-2 weeks). Store at -20 °C 12 months. Store at -80°C long term.
Note	For research use only.

Data Image



SDS PAGE for dsDNA Antibody

SDS-PAGE



Dot blot analysis with dsDNA Antibody (BGT-ANT-35343) at 1µg/ml.

Lane 1: Plasmid (dsDNA)

Second Ab: Goat Anti-Mouse IgG H&L Polyclonal antibody, HRP at 0.1 µg/mL.

Dot Blot