

Escherichia coli O157:H7 FliC/Flagellin Antibody

Summary

Catalog No.	BGT-ANT-35369
Alternative Names	Flagellin, fliC, Z3013, Bacteria, Escherichia coli O157:H7
Host species	Human
Species reactivity	Escherichia coli O157:H7
Form	Liquid
Storage buffer	0.01M PBS, pH 7.4.
Concentration	1 mg/ml
Purity	>95% by SDS-PAGE.
Clonality	Monoclonal
Isotype	IgG1
Applications	ELISA
Target	FliC/Flagellin
Purification	Protein A or G purified from cell culture supernatant.
Endotoxin level	Please contact with the lab for this information.
Accession	Q7DBI0
Stability and Storage	Use a manual defrost freezer and avoid repeated freeze-thaw cycles. Store at +4°C short term (1-2 weeks). Store at -20 °C 12 months. Store at -80°C long term.