

Sulfo-Cy5, NHS Ester

Cat #: BGT-CHM-209

Storage: Store at -20°C protected from light.

Product Introduction

Cy5 NHS Ester are commonly used for the labeling of biomolecules, fluorescence imaging, and other fluorescent bioanalyses. They are widely applied to peptides, proteins, nucleotides, and so on.

Product Properties

Formula: C₃₇H₄₂KN₃O₁₀S₂


MW: 791.99

Solvent: DMSO/DMF/Water

Spectral properties: Ex / Em = 646 / 664 nm

Storage: Store at -20°C away from light and moisture, valid for one year

Structural formula



Note:

The product listed herein is for research use only and is not intended for use in human or clinical diagnosis.