

Human CD103/ITGAE Recombinant Protein

Summary


Catalog No.	BGT-PPT-11081
Alternative Names	Mucosal lymphocyte 1 antigen, Integrin alpha-E, Integrin alpha-IEL, HML-1 antigen, CD103, ITGAE
Form	Lyophilized
Storage buffer	Lyophilized from a solution in PBS pH 7.4, 0.02% Sodium lauroylsarcosinate , 1mM EDTA, 4% Trehalose, 1% Mannitol.
Purity	>90% as determined by SDS-PAGE.
Applications	SDS-PAGE, WB, ELISA, Immunogen, Bioactivity testing in progress
Endotoxin level	Please contact with the lab for this information.
Expression system	E. coli
Accession	P38570
Protein length	Val37-Asn140, N-His- SUMO Tag
Nature	Recombinant
Predicted molecular weight	23.82 kDa
Stability and Storage	Use a manual defrost freezer and avoid repeated freeze thaw cycles.Store at 2 to 8 °C for one week .Store at -20 °C for twelve months from the date of receipt.
Reconstitution	Reconstitute in sterile water for a stock solution.A copy of datasheet will be provided with the products, please refer to it for details.

Species Human

Shipping In general, proteins are provided as lyophilized powder/frozen liquid. They are shipped out with dry ice/blue ice unless customers require otherwise.

Note For research use only.

Data Image



SDS-PAGE

SDS PAGE for Recombinant Human CD103/ITGAE protein