

## Human Calmodulin 2 Recombinant Protein

Cat #: BGT-PPT-34105

Size: 20ug

Storage: Store at -60°C or below. Avoid repeated freeze / thaw cycles.

### Synonyms:

PHKD, CAMII, PHKD2, Calmodulin 2, CALM 2, Calmodulin2, CAM II, LP7057 Protein, PHKD 2, Calmodulin-2, CALM2, Phosphorylase Kinase Delta

### Product Information

**Type:** Recombinant Protein

**Species:** Human

**Source:** Calmodulin 2, 1-149aa, Human , E.coli

**Purity:** > 90% by SDS - PAGE


**Buffer:** Liquid in 20mM Tris pH 7.5

**Sequence:** MADQLTEEQI AEFKEAFSLF DKDGDGTITT KELGTVMRSL GQNPTEAELQ DMINEVDADG NGTIDFPEFL  
TMMARKMKDT DSEEEIREAFRVFDKDGNGY ISAAELRHVM TNLGEKLTDE EVDEMIREAD IDGDGQVNYE EFVQMMTAK

**Molecular Weight:** 16kDa (149aa), confirmed by MALDI-TOF

**Entrez Gene:** 805

**UniProt:** P62158



**Note:**

The product listed herein is for research use only and is not intended for use in human or clinical diagnosis.